



VSR-50

Adjustable angle, dual table capability and up to 100 rocking motions per minute



#517100-00 Dual Stacking Table, doubles capacity

VSR-50[®] Variable Speed Rocker

The VSR-50 Variable Speed Rocker provides a rhythmic side-to-side motion and variable speed control for gentle to vigorous movement. This makes it ideal for gel staining/destaining, Northern and Southern hybridization, Western blotting and membrane rinsing.

The VSR-50 works with either a 12" x 14" (standard) or a 16" x 20" (optional) table. Both tables are easily removable, come equipped with a non-slip neoprene rubber mat, and allow for an angle adjustment up to 15 degrees. For greater capacity without sacrificing counter space, a second 12" x 14" table can be stacked to create a double platform.

The VSR-50's steel housing supports up to 30 lbs while ensuring smooth operation. Features include variable speed control for 5 to 100 rocking motions per minute, a 0-60 minute timer and a "HOLD" feature.

Specifications

Speed	5 - 100 rocking motions per minute (variable)
Electrical	115 VAC, 50/60 Hz; 220 VAC, 50/60 Hz
Tilt angle	Adjustable to 15°
Timer	0 to 60 minutes with "HOLD" feature
Weight	23 lbs
Dimensions	11" x 12" x 5" (d x w x h) without table

Part No. Description

512000-00	VSR-50 Variable Speed Rocker with 12" x 14" table; 115 V, 0.6 A
512500-00	VSR-50 Variable Speed Rocker with 12" x 14" table; 220 V, 0.3 A
517100-00	12" x 14" Stacking Table with 5" standoffs for dual stacking
517200-00	16" x 20" Table

Lab Plus Series Hint

All Lab Plus units are suitable for use in a cold room or incubator. Ambient operating range of 2°C to 40°C

NEXT-DAY DELIVERY ON THIS ITEM and most parts and accessories
 Call 1-800-584-3776 to check availability